

Screening for 34 pesticides in hops using GC-MS/MS

Robert H. Clifford, Nicole Lock, Laura Chambers, William Lipps, Zhuangzhi 'Max' Wang; Shimadzu Scientific Instruments, Columbia, MD, USA

Introduction

From a consumer safety point-of-view, quantitation of the pesticide residues hops has begun to attract wide interest. There are several problems associated with analysis of pesticide residues in hops. First and foremost, there are very few regulatory guidelines established to define which pesticides to include or what the detection limits should be. In fact the US government's 40 CFR Part 180 states individual tolerances must be established for miscellaneous commodities intentionally not included in any group including hops. Secondly the matrix is very complex with significant interferences. Finally, sample load is growing exponentially, so the chosen method must be quick and easy to perform. Trace level pesticide analysis in complex food matrices have been done for many years with similar challenges, thus many of the analytical protocols emerging for the hops matrix are based on these well-established techniques. Triple-quadrupole GC-MS/MS operated in MRM mode provides significant sensitivity and selectivity, but method development can be expensive and time consuming. This poster describes streamlined method development process for analysis of 34 pesticide residues in hops using a QuEChERS sample preparation method, followed by GC-MS/MS detection and quantitation. The pesticides are from 5 classes of compounds including organonitrogen, organophosphorus, organochlorines, carbamates, and synthetic pyrethroids.

A few of the target pesticides were not included in the Smart Pesticide Database. For these compounds, the MRM Optimization Tool was used to automatically determine the optimized MRM transitions and collision energies (CE). Once determined, the new transitions are added to the Smart Pesticide Database. Figure 2 shows the graphic results from the MRM Optimization Tool, with 6 transitions for two of the pesticides.




Figure 2: Optimized Transitions for Two Pesticides Using the MRM Optimization Tool

Experimental

Compound List
For this study 34 pesticides were selected for analysis based on the types of pesticides that are commonly used in hop production. The list includes several different compound classes (Table 1).

Organonitrogen Compounds	Synthetic Pyrethroid Compounds	Organophosphorus Compounds
Bupirimate	Bifenthrin	Chlorpyrifos
Etofenprox	Permethrin	Diazinon
Etridiazole (Terrazole)	Cyfluthrin	Malathion
Fenarimol	Deltamethrin	Mevinphos (Phosdrin)
Flutriafol	Flucythrinate	Phosalone
MGK-264	Lambda-cyhalothrin	Pririmphos methyl
Myclobutanil	Tefluthrin	Carbamates and others
Paclobutrazol	Transfluthrin	Metolaxyl
Penconazole	Organochlorines compounds	2-Phenylphenol
Tebuconazole (Follicur)	Dichlorvos (DDVP)	Vinclozolin
Terbutylazine	Endosulfan sulfate	
Triadimefon	gamma-BHC (Lindane)	
Triadimenol (Baytan)	p,p'-DDT	

Table 1: Selected Pesticide Compound Classes Included Organonitrogens, Synthetic Pyrethroids, Organochlorines, Organophosphates, and Carbamates

Method Development

The most difficult part of any triple quadrupole method development process, is determination and optimization of the Multiple Reaction Monitoring (MRM) transitions and collision energies (CE). For this study, the Shimadzu Smart Pesticide Database was used as the foundation for creating the MRM analysis method. The Smart Pesticide Database includes up to six fully optimized MRM transitions and CEs for 479 pesticides and Retention Indices (RI) for accurately predicting compound retention times. The transitions and CEs in the database were optimized using the Shimadzu GCMS-TQ8040 triple quadrupole GC-MS/MS. Figure 1 shows a portion of the Smart Pesticide Database and the method, compound, and transition information.

Serial	Type	Acq. Method	Method No.	Compound Name (S)	Ret. Index (I)	CE (eV)							
1	Target	MSM	1	Atrazine	887	0	0	0	0	0	0	0	0
2	Target	MSM	1	Chlorpyrifos	1027	108	0	0	0	0	0	0	0
3	Target	MSM	1	Alachlor	1154	0	0	0	0	0	0	0	0
4	Target	MSM	1	Chlorpyrifos	1164	0	0	0	0	0	0	0	0
5	Target	MSM	1	Chlorpyrifos	1201	1027	44	0	0	0	0	0	0
6	Target	MSM	1	Malathion	1240	1026	32	0	0	0	0	0	0
7	Target	MSM	1	Chlorpyrifos	1248	1027	44	0	0	0	0	0	0
8	Target	MSM	1	Malathion	1280	1026	32	0	0	0	0	0	0
9	Target	MSM	1	Endosulfan	1296	89	11	0	0	0	0	0	0
10	Target	MSM	1	Dichlorvos	1348	1194	66	0	0	0	0	0	0
11	Target	MSM	1	Chlorpyrifos	1398	1027	44	0	0	0	0	0	0
12	Target	MSM	1	Malathion	1400	1026	32	0	0	0	0	0	0

Figure 1: Example of Information Found in the Smart Pesticide Database Used to Create an MRM Analysis Method

Serial	Compound Name	Ret. Index (I)	Ret. Time	Type	Ion1	CE	Ratio	Type	Ion2	CE	Ratio	Type	Ion3	CE	Ratio
1	Dichlorvos	1252	4.565	T	109.05-79.10	7	100.00	Ref1	165.05-63.10	13	44.15	Ref2	210.95-163.00	5	10.19
2	Mevinphos	1427	5.642	T	127.05-109.00	11	100.00	Ref1	192.05-127.00	13	47.34	Ref2	127.05-99.00	15	35.24
3	Endosulfan	1459	5.891	T	210.95-140.00	15	96.56	Ref1	182.55-140.00	15	96.56	Ref2	210.95-140.00	15	96.56
4	2-Phenylphenol	1533	6.483	T	109.10-141.10	13	100.00	Ref1	109.10-115.10	25	91.99	Ref2	170.10-141.10	23	86.19
5	Lindane	1779	8.601	T	180.15-148.10	15	99.99	Ref1	180.15-148.10	15	99.99	Ref2	210.95-163.00	19	33.81
6	Terbutylazine	1782	8.604	T	229.10-173.10	7	100.00	Ref1	214.10-171.10	19	78.34	Ref2	214.10-132.10	9	59.16
7	Diazinon	1790	8.766	T	104.10-179.20	13	100.00	Ref1	248.05-152.10	7	61.75	Ref2	248.05-137.10	17	61.34
8	Tefluthrin	1816	9.002	T	177.05-127.10	17	100.00	Ref1	177.05-127.10	17	100.00	Ref2	197.05-141.10	13	31.28
9	Malathion	1894	9.730	T	217.05-157.00	15	100.00	Ref1	217.05-157.00	15	100.00	Ref2	280.05-212.00	15	71.11
10	Transfluthrin	1903	9.815	T	163.05-127.10	7	100.00	Ref1	163.05-91.10	15	82.75	Ref2	163.05-143.00	17	75.90
11	Metolaxyl	1915	9.926	T	234.10-148.10	19	100.00	Ref1	234.10-174.10	11	75.22	Ref2	209.15-190.20	9	64.50
12	Phosphos methyl	1941	10.167	T	200.10-125.10	23	100.00	Ref1	200.10-203.10	11	53.89	Ref2	230.05-223.00	17	54.23
13	Malathion	1964	10.377	T	127.05-99.10	7	100.00	Ref1	173.10-127.10	13	66.84	Ref2	173.10-127.10	7	64.75
14	Chlorpyrifos	1980	10.529	T	113.95-257.80	19	100.00	Ref1	215.95-259.90	19	74.59	Ref2	285.95-257.90	9	47.29
15	Triadimefon	2003	10.738	T	208.05-111.10	23	100.00	Ref1	208.05-127.10	15	89.54	Ref2	210.95-183.10	9	43.88
16	MGK-264	2008	10.890	T	164.10-99.10	13	100.00	Ref1	164.10-99.10	13	100.00	Ref2	164.10-99.10	28	55.15
17	Penconazole	2063	11.283	T	248.10-157.10	25	100.00	Ref1	159.05-123.10	19	50.14	Ref2	248.10-192.10	15	45.77
18	Triadimefon	2092	11.541	T	168.15-70.00	9	100.00	Ref1	128.05-65.10	23	38.42	Ref2	112.05-58.10	11	27.68
19	Paclobutrazol	2124	11.899	T	226.05-125.10	11	100.00	Ref1	226.05-167.10	9	31.19	Ref2	248.05-127.10	11	32.42
20	Phosdrin	2125	12.104	T	123.05-95.10	13	100.00	Ref1	119.05-123.10	15	65.69	Ref2	123.05-75.10	25	52.38
21	Myclobutanil	2200	12.502	T	179.05-125.00	15	100.00	Ref1	179.05-152.00	9	55.34	Ref2	179.05-96.10	29	56.16
22	Bupirimate	2204	12.535	T	273.10-108.10	15	100.00	Ref1	273.10-193.10	7	67.72	Ref2	193.15-81.10	25	74.80
23	Endosulfan sulfate	2269	13.865	T	271.80-228.80	21	100.00	Ref1	271.80-224.80	17	22.29	Ref2	271.80-141.00	11	22.11
24	p,p'-DDT	2367	15.919	T	235.05-165.20	29	100.00	Ref1	237.05-165.20	23	64.85	Ref2	235.05-199.10	17	13.44
25	Tebuconazole	2399	14.184	T	250.10-125.10	19	100.00	Ref1	250.10-70.10	9	40.63	Ref2	252.10-127.10	23	55.38
26	Bifenthrin	2471	14.707	T	181.15-166.10	13	100.00	Ref1	181.15-166.10	27	90.00	Ref2	166.10-164.20	29	4.99
27	Phosalone	2566	16.432	T	162.05-111.00	15	100.00	Ref1	162.05-75.10	27	53.27	Ref2	162.05-184.00	9	38.67
28	lambda-Cyhalothrin	2597	15.748	T	197.05-141.10	11	100.00	Ref1	208.10-181.10	7	97.01	Ref2	197.05-161.10	7	54.32
29	Penconazole	2631	16.601	T	251.05-139.00	15	100.00	Ref1	251.05-111.10	29	42.14	Ref2	330.05-139.10	9	34.45
30	Phenothiazin	2709	16.562	T	183.10-163.10	15	100.00	Ref1	183.10-168.10	15	107.11	Ref2	163.05-127.10	7	106.13
31	Cyfluthrin	2793	17.202	T	226.05-206.10	13	100.00	Ref1	199.10-170.10	25	70.95	Ref2	206.05-151.10	19	64.85
32	Bifenthrin	2870	17.812	T	163.15-158.10	11	100.00	Ref1	163.15-167.10	17	89.29	Ref2	176.20-163.20	11	5.78
33	Phenylazine	2876	17.860	T	199.10-157.10	9	100.00	Ref1	199.10-167.10	23	94.17	Ref2	225.10-119.10	19	18.37
34	Deltamethrin	2961	19.900	T	252.05-93.10	19	100.00	Ref1	161.10-152.10	23	87.80	Ref2	252.05-172.00	7	50.01

GC-2010 Plus	GCMS-TQ8040
Injection Port	250 °C
Column	1 µL splitless injection, 1 minute sampling time
Oven Program	SH-Rxi-5Sil MS, 30 m x 0.25 mm x 0.25 µm film Helium carrier gas Constant Linear Velocity mode, 40.0 cm/second 85 °C (hold 1 minute) 25 °C/minute to 160 °C 10 °C/minute to 240 °C 10 °C/minute to 290 °C (hold 3 minutes)
Transfer Line	300 °C
Mass Spectrometer	GCMS-TQ8040
Acquisition Mode	MRM
Ion Source	230 °C Electron ionization mode, 70 eV
Collision Gas	Argon, 200 kPa
MRM Loop Time	Optimized with Smart MRM

Table 3 Optimized Instrument Conditions for Analysis of Pesticides in Hop Samples using the Shimadzu GCMS-TQ8040

Sample Preparation - QuEChERS	QuEChERS Extraction Steps Followed by Cartridge SPE Cleanup
2 g ground hops into 50 mL centrifuge tube	Add 10 mL water
	Add 10 mL acetonitrile
	Shake/vortex 30 minutes
	Add QUENCHERS salt
	Cap and shake 1 minute
	Centrifuge 5 minutes
	Add 10.5 mL MgCl ₂ 10.6 mL SPE cartridge into 50 mL vial
	Shake with 20 mL acetonitrile
	Load 1 mL QuEChERS extract onto cartridge
Flute with 15 mL 1:1 acetonitrile	Concentrate eluate to 4-6 µL




Figure 3: The MRM Analysis Method is Created Automatically and Optimized for Sensitivity

MRM Optimization Tool

Figure 4: Workflow using the MRM Optimization Tool

Information used to create the analysis method is shown in Table 2. It includes a compound table, retention indices and retention times, one transition with optimized CE for quantitation, and two reference transitions. Area ratios are also empirically determined, and can be used as part of the laboratory QA/QC program.




Figure 5 Exponential Calibration Curves for Two Pesticides, 1 to 200 ppb




Figure 6 Example of Overlaid MRM Chromatograms For Two Pesticides in the 1-ppb Calibration Standard

Sample Repeatability

Two different hops samples were processed using the QuEChERS procedure. The extracts were spiked with the pesticide mix at 25 ppb and analyzed in triplicate using the optimized MRM method. Chromatograms in Figure 7 illustrate how the MRM technique can be used to select the target compound from a complex matrix background, and produce reliable, reproducible results at low concentrations.




Figure 7 MRM Chromatograms of Two Hops Samples Spiked at 25 ppb and Analyzed in Triplicate

Summary and Conclusion

The data presented illustrate how a triple quadrupole GC-MS/MS operated in the MRM mode, can be used to analyze for trace-level pesticide residues in complex plant matrices such as hops. The matrix was extracted using a QuEChERS kit, and interferences removed using an SPE cartridge. The resulting extracts were analyzed in triplicate using MRM transitions provided in the Smart Pesticide Database or individually optimized using the MRM Optimization Tool, with repeatability of 6% or better. The MRM method was fully optimized in just a few minutes, target compounds were selectively identified against the co-eluting matrix interferences, and quantitated at the parts-per-billion range.

Acknowledgement

The authors wish to acknowledge Restek Corporation, Bellefonte, PA for useful discussions and advice regarding column selection, standards, and providing the QuEChERS materials used in the development of this method.