

HPE-9

*American Society of Brewing Chemists
HOPS CHECK SERVICE*

9/30/2015

Subscriber Number	<i>Spectrophotometric Methods</i>			<i>Conductometric Methods</i>			<i>HPLC Method</i>				Standard
	<i>Alpha Acids (%)</i>	<i>Beta Acids (%)</i>	<i>Hop Storage Index</i>	<i>ASBC C.V. (%)</i>	<i>EBC LCV (%)</i>	<i>MEBAK LCV (%)</i>	<i>Alpha Acids (%)</i>	<i>Beta Acids (%)</i>	<i>Cohumulone (% of total)</i>	<i>Colupulone (% of total)</i>	
B 5	66.70	17.72	0.235			63.80	60.26	16.39	28.56	52.290	ICE-3
B 6	65.37	15.47	0.251				58.42	15.98	28.49	52.800	
B 8	64.71	16.87	0.234			61.28	59.38	16.18	28.53	53.010	ICE-3
B 10	64.76	17.62	0.229				58.20	16.87	28.00	51.230	ICE-3
B 17	68.48	15.50	0.260								
B 18							59.68	16.27	* 3.67	* 28.200	ICE-3
B 115	66.46	17.05	0.241			62.21	59.84	16.29	28.53	52.670	ICE-3
B 127	64.61	15.18	0.242			62.34	59.66	15.55	* 26.68	51.450	ICE-3
B 142	67.07	16.72	0.255				59.04	15.86	28.98	52.620	ICE-3
B 178							* 54.40	15.10	28.50	52.600	ICE-3
B 196	65.75	15.34	0.270								
N	9	9	9			4	8	9	7	8	
Mean	65.99	16.39	0.246			62.41	59.31	16.05	28.51	52.334	
Min	64.61	15.18	0.229			61.28	58.20	15.10	28.00	51.230	
Max	68.48	17.72	0.270			63.80	60.26	16.87	28.98	53.010	
Std	1.30	1.02	0.014			1.04	0.71	0.51	0.28	0.648	

* Not included in mean.

Confidentiality Notice: The data contained herein are confidential and intended for subscriber use only. Data from other subscribers are included for comparison purposes only.