

HPW-10

American Society of Brewing Chemists
HOPS CHECK SERVICE

10/31/2015

Subscriber Number	<i>Spectrophotometric Methods</i>			<i>Conductometric Methods</i>			<i>HPLC Method</i>				Standard
	<i>Alpha Acids (%)</i>	<i>Beta Acids (%)</i>	<i>Hop Storage Index</i>	<i>ASBC C.V. (%)</i>	<i>EBC LCV (%)</i>	<i>MEBAK LCV (%)</i>	<i>Alpha Acids (%)</i>	<i>Beta Acids (%)</i>	<i>Cohumulone (% of total)</i>	<i>Colupulone (% of total)</i>	
B 6	6.58	4.33	0.438				6.25	4.08	25.83	49.660	
B 8	6.40	4.64	0.426		6.48	7.22	6.35	4.10	25.29	50.680	ICE-3
B 10	6.32	4.89	0.434			6.95	6.37	4.40	26.33	47.930	ICE-3
B 17	6.53	4.02	0.482								
B 99	6.36	4.83	0.450		6.77	7.70					
B 115	6.25	4.58	0.467			7.60	6.52	4.10	25.70	50.000	ICE-3
B 127	6.02	4.30	0.480								
B 142	6.62	4.69	0.491				6.01	3.77	26.67	49.370	ICE-3
B 169	5.79	4.73	* 0.550				5.83	3.88	26.15	47.970	ICE-3
B 171	6.57	4.40	0.445								
B 178	5.97	3.99	0.419				6.00	3.80	25.80	50.900	ICE-3
B 196	6.46	4.64	0.450								
B 200	6.20	4.60	0.472								
N	13	13	12		2	4	7	7	7	7	
Mean	6.31	4.51	0.454		6.63	7.37	6.19	4.02	25.97	49.501	
Min	5.79	3.99	0.419		6.48	6.95	5.83	3.77	25.29	47.930	
Max	6.62	4.89	0.491		6.77	7.70	6.52	4.40	26.67	50.900	
Std	0.26	0.28	0.024		0.21	0.35	0.25	0.22	0.45	1.187	

* Not included in mean.

Confidentiality Notice: The data contained herein are confidential and intended for subscriber use only. Data from other subscribers are included for comparison purposes only.